

Infectious diseases *Enfermedades infecciosas*

GDH, 2 A-B diff

Immunochromatographic test for the diagnosis of Infectious diseases in stool samples
Test inmunocromatográfico para el diagnóstico de enfermedades infecciosas en heces



www.operon.es - sales@operon.es - Tel.: +34 976 503 597 - Spain



Clostridium difficile




Name Nombre	Analyte Analito	Sample Muestra	Sensitivity Sensibilidad	Specificity Especificidad	Stick format	Simple format
GDH	Clostridium difficile	Stool samples Heces	96%	98%	X	X
2 A-B diff	Clostridium difficile	Stool samples Heces	94,9%	96,1%	X	X


Sample preparation / Preparación de muestra



Stick GDH




Simple GDH



Stick 2 A-B diff

INDEPENDENT bands for toxin A and toxin B



Simple 2 A-B diff

