



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	108C.USE
GENERAL DESCRIPTION	TRANSYSTEM AMIES A/PL 500PKG
DETAILED DESCRIPTION	1 box containing 10 Vi-pak each containing 50 film-film pouches with 1 tube with 5 ml AMIES Medium with white cap and 1 plastic applicator with blue cap

PRODUCT COMPONENTS

MEDIUM		COMPOSITION
AMIES WITHOUT CHARCOAL		Sodium chloride Potassium Chloride Calcium Chloride Monopotassium Phosphate Magnesium chloride Disodium Phosphate Sodium Thioglycollate Bacto Agar Distilled Water

APPLICATOR	PLASTIC APPLICATOR	SPECIFICATIONS
		TOTAL LENGTH: 147.00 mm $\pm$ 1.00 mm STICK DIAMETER: 2.50 mm $\pm$ 0.20 mm FLOCKED LENGTH: 15.00 mm $\pm$ 2.00 mm FLOCKED DIAMETER: 4.50 mm $\pm$ 0.50 mm STICK MATERIAL: POLYSTYRENE

CAP	BLUE CAP	SPECIFICATIONS
		EXTERNAL CAP DIAMETER: 12.00 mm $\pm$ 0.10 mm CAP LENGTH: 32.80 mm $\pm$ 0.20 mm INSERTION DIAMETER: 9.30 mm $\pm$ 0.05 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.02-07/05/2019) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

TRANSPORT TUBE	TUBE	SPECIFICATIONS
		TUBE LENGTH: 154.50 mm $\pm$ 1.00 mm INSERTION DIAMETER: 9.00 mm $\pm$ 0.50 mm EXTERNAL DIAMETER: 12.00 mm $\pm$ 0.50 mm MATERIAL: POLYPROPYLENE
CAP	WHITE CAP W/O RECEPTACLE	SPECIFICATIONS
		EXTERNAL DIAMETER: 12.00 mm $\pm$ 0.10 mm CAP LENGTH: 32.80 mm $\pm$ 0.20 mm INSERTION DIAMETER: 9.00 mm $\pm$ 0.05 mm MATERIAL: HIGH DENSITY POLYETHYLENE

PACKAGING SPECIFICATION

OUTER BOX QTY	500
PACKAGING	FILM-FILM POUCH
VI-PAK QUANTITY	50

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	20 MONTHS
MEDIUM pH	7.2 $\pm$ 0.3

This bulletin (Rev.02-07/05/2019) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.