



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	147C
GENERAL DESCRIPTION	VIRO-SWAB T/LIQ.PLAS.
DETAILED DESCRIPTION	10 boxes contain 25 film-film pouches each contains 1 plastic applicator with polyester tip and pink cap, and a tube with white cap without receptacle and sponge with 1.2 ml of medium for viral infections (Chlamydia)

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
VIRAL MEDIUM	BUFFERING AGENTS HBSS OSMOTIZED WATER

APPLICATOR	PLASTIC APPLICATOR	SPECIFICATIONS
		TOTAL LENGTH: 147.00 mm ± 2.00 mm STICK DIAMETER: 2.50 mm ± 0.10 mm FLOCKED LENGTH: 16.00 mm ± 3.00 mm FLOCKED DIAMETER: 4.50 mm ± 0.50 mm STICK MATERIAL: POLYSTYRENE
TRANSPORT TUBE	TUBE	SPECIFICATIONS
		TUBE LENGTH: 154.50 mm ± 1.00 mm INSERTION DIAMETER: 9.00 mm ± 0.50 mm EXTERNAL DIAMETER: 12.00 mm ± 0.50 mm MATERIAL: POLYPROPYLENE

This bulletin (Rev.07-02/11/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

SPONGE	SPONGE	SPECIFICATIONS
		LENGTH: 37.50 mm + 2.00 mm -1.50 mm WIDTH: 10.00 mm ± 0.50 mm THICKNESS: 5.00 mm ± 0.50 mm MATERIAL: POLYURETHANE
CAP	WHITE CAP W/O RECEPTACLE	SPECIFICATIONS
		EXTERNAL DIAMETER: 12.00 mm ± 0.10 mm CAP LENGTH: 32.80 mm ± 0.20 mm INSERTION DIAMETER: 9.00 mm +0.05 mm MATERIAL: HIGH DENSITY POLYETHYLENE
CAP	PINK CAP	SPECIFICATIONS
		EXTERNAL CAP DIAMETER: 12.00 mm ± 0.20 mm CAP LENGTH: 32.90 mm ± 0.50 mm INSERTION DIAMETER: 9.30 mm + 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

PACKAGING SPECIFICATION

INNER BOX QTY	25
OUTER BOX QTY	250
PACKAGING	FILM-FILM POUCH

This bulletin (Rev.07-02/11/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	15 MONTHS
MEDIUM pH	7.2 ± 0.3