



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	150C
GENERAL DESCRIPTION	TAMPONE SECCO A/LEGNO 1000PKG
DETAILED DESCRIPTION	1 outer box contain 10 inner boxes each composed of 100 wooden swabs in tube with red cap.

PRODUCT COMPONENTS

APPLICATOR	WOODEN APPLICATOR	SPECIFICATIONS
		TOTAL LENGTH: 150.50 mm \pm 2.00 mm STICK DIAMETER: 2.30 mm FIBER TIP DIAMETER: 4.80 mm \pm 1.00 mm FIBER TIP LENGTH: 15.00 mm \pm 2.00 mm STICK MATERIAL: WOOD
CAP	RED CAP	SPECIFICATIONS
		EXTERNAL CAP DIAMETER: 12.00 mm \pm 0.20 mm CAP LENGTH: 32.90 mm \pm 0.50 mm INSERTION DIAMETER: 9.30 mm + 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE



TECHNICAL DATA SHEET

TRANSPORT TUBE	TUBE	SPECIFICATIONS
		TUBE LENGTH: 154.50 mm \pm 1.00 mm INSERTION DIAMETER: 9.00 mm \pm 0.50 mm EXTERNAL DIAMETER: 12.00 mm \pm 0.50 mm MATERIAL: POLYPROPYLENE

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	BULK

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS