



## TECHNICAL DATA SHEET

### PRODUCT DESCRIPTION

<b>PRODUCT CODE</b>	155C
<b>GENERAL DESCRIPTION</b>	TAMPONE SECCO A/PL 1000 PKG
<b>DETAILED DESCRIPTION</b>	1 outer box contain 10 inner boxes each composed of 100 fibre tip swabs in tube with white cap.

### PRODUCT COMPONENTS

<b>APPLICATOR</b>	PLASTIC APPLICATOR	<b>SPECIFICATIONS</b>
		STICK DIAMETER: 2.50 mm $\pm$ 0.20 mm FIBRE TIP LENGTH: 15.00 mm $\pm$ 2.00 mm FIBRE TIP DIAMETER: 4.50 mm $\pm$ 0.50 mm TOTAL LENGTH: 147.00 mm $\pm$ 2.00 mm MATERIAL: POLYSTYRENE
<b>CAP</b>	WHITE CAP	<b>SPECIFICATIONS</b>
		EXTERNAL CAP DIAMETER: 12.00 mm $\pm$ 0.10 mm CAP LENGTH: 32.80 mm $\pm$ 0.20 mm INSERTION DIAMETER: 9.30 mm + 0.05 mm MATERIAL: HIGH DENSITY POLYETHYLENE



## TECHNICAL DATA SHEET

TRANSPORT TUBE	TUBE	SPECIFICATIONS
		TUBE LENGTH: 154.50 mm $\pm$ 1.00 mm INSERTION DIAMETER: 9.00 mm $\pm$ 0.50 mm EXTERNAL DIAMETER: 12.00 mm $\pm$ 0.50 mm MATERIAL: POLYPROPYLENE

### **PACKAGING SPECIFICATION**

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	BULK

### **IMPORTANT INFO**

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS