



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	305C
GENERAL DESCRIPTION	UTM 3ML RT PERN.FLOC AP 300PKG
DETAILED DESCRIPTION	1 outer box contain 6 inner boxes each composed of 50 paper-film pouches each containing 1 flocked swab individually packaged in paper-film pouch and 1 tube with red cap filled with 3mL of UTM solution.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	CONICAL TUBE	SPECIFICATIONS
		TUBE LENGTH: 101.50 mm ± 0.20 mm TUBE DIAMETER: 16.30 mm ± 0.20 mm INSIDE DIAMETER: 14.10 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		LOWER CAP DIAMETER: 21.50 mm ± 0.20 mm CAP LENGTH: 16.30 mm ± 0.20 mm HIGHER CAP DIAMETER: 20.50 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.09-29/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	FLOCKED SWAB	SPECIFICATIONS
		BREAK POINT: 98.50 mm \pm 1.00 mm STICK DIAMETER: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 152.00 mm \pm 2.00 mm FLOCKED DIAMETER: 4.00 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 \pm 0.2

This bulletin (Rev.09-29/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.