



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	346C
GENERAL DESCRIPTION	UTM RT 3ML+REG/FLOCK PF100MM
DETAILED DESCRIPTION	6 boxes containing 50 paper film pouches each containing 1 tube with red cap and 3 ml of UTM and 1 regular flocked swab paper film pouched.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	CONICAL TUBE	SPECIFICATIONS
		TUBE LENGTH: 101.50 mm ± 0.20 mm TUBE DIAMETER: 16.30 mm ± 0.20 mm INSIDE DIAMETER: 14.10 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		LOWER CAP DIAMETER: 21.50 mm ± 0.20 mm CAP LENGTH: 16.30 mm ± 0.20 mm HIGHER CAP DIAMETER: 20.50 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.07-16/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	FLOCKED REGULAR MARKED SWAB	SPECIFICATIONS
		STICK DIAMETER: 2.50 mm \pm 0.10 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm FLOCKED DIAMETER: 5.50 mm \pm 0.50 mm TOTAL LENGTH: 151.00 mm \pm 2.00 mm BREAK POINT: 98.50 mm \pm 1.00 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 \pm 0.2

This bulletin (Rev.07-16/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.