



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	348C
GENERAL DESCRIPTION	UTM 2,5ML BULK NO BEADS
DETAILED DESCRIPTION	1 outer box containing 6 inner boxes each comped of 50 tubes with red cap filled UTM Medium

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	CONICAL TUBE	SPECIFICATIONS
		TUBE LENGTH: 101.50 mm ± 0.20 mm TUBE DIAMETER: 16.30 mm ± 0.20 mm INSIDE DIEMETER: 14.10 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	CAP RED S \ RECEPTACLE	SPECIFICATIONS
		HEIGHT CAP: 16.30 mm ± 0.20 mm CAP DIAMETER: 21.50 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.01-17/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	BULK

IMPORTANT INFO

STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7,30 ± 0,20