



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	361C
GENERAL DESCRIPTION	UTM MINI 1ML RT+MINITIP FLOC.
DETAILED DESCRIPTION	6 boxes containing 50 paper-film pouches each containing 1 tube with red screw cap and 1 ml of UTM and 1 urethral flocked swab paper film pouched.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm ± 0.20 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.06-16/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	FLOCKED SWAB	SPECIFICATIONS
		STICK DIAMETER: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 151.00 mm \pm 2.00 mm FLOCKED DIAMETER: 3.50 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm BREAK POINT: 78.50 mm \pm 1.00 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 \pm 0.2