



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	502CS01
GENERAL DESCRIPTION	REGULAR APPL. NY/VEL B/01PCS
DETAILED DESCRIPTION	1 outer box contain 10 inner boxes each composed of 100 regular flocked swabs individually wrapped in paper-film pouch.

PRODUCT COMPONENTS

		SPECIFICATIONS
APPLICATOR	REGULAR FLOCKED SWAB	TOTAL LENGTH: 151.00 mm ± 2.00 mm STICK DIAMETER: 2.50 mm ± 0.10 mm BREAKING POINT: 78.50 mm ± 1.00 mm FLOCKED LENGTH: 16.00 mm ± 2.00 mm FLOCKED DIAMETER: 5.50 mm ± 0.50 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	ETO TREATMENT
STORAGE TEMPERATURE	2/30°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.09-16/11/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.