



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	551C
GENERAL DESCRIPTION	MINITIP PL/APP NY/VEL DRY TUBE
DETAILED DESCRIPTION	1 outer box contain 10 inner boxes each composed of 100 flocked swabs in tube with green cap.

PRODUCT COMPONENTS

APPLICATOR	FLOCKED SWAB	SPECIFICATIONS
		STICK DIAMETER: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 151.00 mm \pm 2.00 mm FLOCKED DIAMETER: 3.50 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm BREAK POINT: 78.50 mm \pm 1.00 mm STICK MATERIAL: POLYSTYRENE
TRANSPORT TUBE	TUBE	SPECIFICATIONS
		TUBE LENGTH: 155.00 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.5 mm EXTERNAL TUBE DIAMETER: 12.0 mm \pm 0.5 mm MATERIAL: POLYPROPYLENE



TECHNICAL DATA SHEET

CAP	GREEN CAP	SPECIFICATIONS
		CAP LENGTH: 32.80 mm \pm 0.50 mm INSERTION DIAMETER: 9.30 mm + 0.20 mm EXTERNAL CAP DIAMETER: 12.00 mm \pm 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE
CAP + TUBE	CAP + TUBE	SPECIFICATIONS
		TOTAL LENGTH: 175.50 mm \pm 3.00 mm

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	BULK

IMPORTANT INFO

STERILIZATION METHOD	ETO TREATMENT
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.05-28/02/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.