



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

| | |
|-----------------------------|--|
| PRODUCT CODE | 553C |
| GENERAL DESCRIPTION | FLEXIBLE FLOQSWAB IN DRY TUBE BP100 |
| DETAILED DESCRIPTION | 1 outer box contain 10 inner boxes each composed of 100 flocked swabs in tube with orange cap. |

PRODUCT COMPONENTS

| | | |
|-----------------------|--------------|---|
| APPLICATOR | FLOCKED SWAB | SPECIFICATIONS |
| | | BREAK POINT: 98.50 mm \pm 1.00 mm STICK DIAMETER: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 152.00 mm \pm 2.00 mm FLOCKED DIAMETER: 4.00 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm STICK MATERIAL: POLYSTYRENE |
| TRANSPORT TUBE | TUBE | SPECIFICATIONS |
| | | TUBE LENGTH: 155.00 mm \pm 1.00 mm INSERTION DIAMETER: 9.0 mm \pm 0.5 mm EXTERNAL TUBE DIAMETER: 12.0 mm \pm 0.5 mm MATERIAL: POLYPROPYLENE |

This bulletin (Rev.07-28/02/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

| | | |
|-------------------|------------|---|
| CAP | ORANGE CAP | SPECIFICATIONS |
| | | EXTERNAL CAP DIAMETER: 12.00 mm \pm 0.20 mm CAP LENGTH: 32.90 mm \pm 0.50 mm INSERTION DIAMETER: 9.30 mm + 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE |
| CAP + TUBE | CAP + TUBE | SPECIFICATIONS |
| | | TOTAL LENGTH: 175.50 mm \pm 3.00 mm |

PACKAGING SPECIFICATION

| | |
|----------------------|------|
| INNER BOX QTY | 100 |
| OUTER BOX QTY | 1000 |
| PACKAGING | BULK |

IMPORTANT INFO

| | |
|-----------------------------|------------------------|
| STERILIZATION METHOD | ETO TREATMENT |
| STORAGE TEMPERATURE | 2/30°C – DO NOT FREEZE |
| SHELF LIFE | 36 MONTHS |

This bulletin (Rev.07-28/02/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.